

# DNC SYSTEM SIMPLE INSTRUCTION

To Memory or DNC from top

DNC from mid

{ % O7790 (READ) or (USBREAD) (NAME) (1234.NC) Your File name % <div style="border: 1px solid black; border-radius: 15px; padding: 5px; width: fit-content; margin: 10px auto;"><p>Alternate command. /T0001 as (READ) /T0021 as (USBREAD) /T0003 as (NAME)</p></div>	{ % O7791 (READ) or (USBREAD) (NAME) (1234.NC) Your File name (NNO) (100) Sequence Number % }
--	--

☞ Pulling Data { PUNCH O7790  
then  
READ

☞ Saving Data

Simply punch the program.

-> If "O1234", "O1234.NCT" will be created.